

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Human Cardiac Troponin-I Monoclonal Antibody***

Catalog No.	10030
Description	Purified Human Cardiac Troponin-I Monoclonal Antibody
Purity	>95% by SDS-PAGE. Purified by protein A agarose chromatography
Protein	1 – 10 mg/ml, dependent on production batch
Buffer	PBS, pH 7.4, containing 0.1% sodium azide as preservative
Storage	2-8°C short term, -20 °C or colder long term
Immunogen:	Highly purified human cardiac Troponin-I
Cross reactivity	This antibody reacts with human cardiac Troponin I. There is no cross reactivity human skeletal Troponin I

P.O. Box 5205 • West Chester, PA 19380 • USA