

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Goat Anti-Human Cardiac Troponin-I Peptide-4 Antibody

Catalog No.	10050
Description	Goat anti-TNI peptide 4, affinity purified (-RGEKGRALSTRCQPLELA, aa 69-86)
Purity	>95% pure. Affinity purified with solid phase human cardiac troponin I peptide 4
Protein	1-10 mg/ml (E1% at 280 nm = 13.0)
Buffer	10 mM Sodium Phosphate, 0.85% (w/v) NaCl, pH 7.2, 0.05% NaN ₃
Storage	2-8°C
Immunogen	Highly purified human cardiac Troponin-I
Comments	Cross reactivity with human skeletal Troponin I is less than 1% by ELISA.

P.O. Box 5205 • West Chester, PA 19380 • USA