

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Human Cardiac FABP Monoclonal Antibody

Catalog No.	10090
Description	Purified human cardiac fatty acid binding protein monoclonal antibody (clone 31A15.1)
Subclass	IgG 2b, kappa
Purity	>95% by SDS-PAGE. Purified by protein A agarose chromatography
Protein	1-10 mg/ml
Buffer	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2, containing 0.05% sodium azide as preservative
Storage	2-8°C
Immunogen:	Purified human cardiac fatty acid binding protein

P.O. Box 5205 • West Chester, PA 19380 • USA