

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818 • [Info@lifediagnostics.com](mailto:Info@lifediagnostics.com)  
[www.lifediagnostics.com](http://www.lifediagnostics.com)

## PRODUCT DATA SHEET

**PRODUCT:** **Mouse Serum Control**

**CATALOG #:** MSCON01

**LOT #:** N-A2209D

**FORMAT:** Lyophilized mouse serum. Reconstitute by adding 200 µl of distilled or deionized water.

**DESCRIPTION:** Contains the indicated levels of serum proteins (see table below)\*

**STORAGE:** Store lyophilized control at -20°C or colder

<b>Analyte</b>	<b>Concentration</b>
Alpha-Macroglobulin	3.37 +/- 0.92 mg/ml
C-Reactive Protein	116 +/- 29 µg/ml
Clusterin	12.8 +/- 1.7 µg/ml
Haptoglobin	238 +/- 12 µg/ml
Hemopexin	1.63 +/- 0.41 mg/ml
IgG	1.96 +/- 0.15 mg/ml
IgM	0.87 +/- 0.21 mg/ml