

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Goat Anti Cat $\alpha$ -1-Acid Glycoprotein (Affinity Purified)***

Catalog No.	18132
Description	Goat anti-cat $\alpha$ -1-acid glycoprotein, affinity pure IgG
Purity	>90%, affinity purified on cat $\alpha$ -1-acid glycoprotein agarose
Protein	1-10 mg/ml (Assuming $A_{280} = 14$ for a 10 mg/ml solution)
Buffer	10 mM sodium phosphate, 150 mM NaCl, pH 7.2, 0.1 % $\text{NaN}_3$
Immunogen	Pure cat $\alpha$ -1-acid glycoprotein
Applications	ELISA (as capture and/or detection antibody), Western blot
Storage	2-8°C short term, -20°C long term

P.O. Box 5205 • West Chester PA 19380 • USA