

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Rabbit Anti Mouse α -1-Acid Glycoprotein (Affinity Purified)

Catalog No.	18140
Description	Rabbit anti-mouse α -1-acid glycoprotein
Purification	Affinity purified on rat α -1 acid glycoprotein agarose
Applications	ELISA, Western blot
Protein	0.5 – 2.5 mg/ml (Assuming $A_{280} = 14$ for a 10 mg/ml solution)
Buffer	10 mM sodium phosphate, 150 mM NaCl, pH 7.2, 0.05 % NaN_3
Storage	2-8°C
Source	Serum from rabbits immunized with purified mouse α -1-acid glycoprotein

P.O. Box 5205 • West Chester, PA 19380 • USA