

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Rabbit Anti Mouse Haptoglobin (Affinity Purified)

| | |
|--------------|--|
| Catalog No. | 18141 |
| Description | Rabbit anti-mouse haptoglobin |
| Purity | Affinity purified on mouse haptoglobin agarose |
| Applications | ELISA , Western blot |
| Protein | 0.5 – 5 mg/ml (Assuming $A_{280} = 14$ for a 10 mg/ml solution) |
| Buffer | 10 mM sodium phosphate, 150 mM NaCl, pH 7.2, 0.1 % NaN_3 |
| Storage | 2-8°C |
| Source | Serum from rabbits immunized with purified mouse haptoglobin |

P.O. Box 5205 • West Chester, PA 19380 • USA