

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Rabbit Anti Rat Alpha-2-Macroglobulin (Affinity Purified)***

Catalog No.	18142
Purification	Affinity purified on rat alpha-2-macroglobulin agarose
Applications	ELISA, Western blot
Protein	0.25 – 2.5 mg/ml (Assuming $A_{280} = 14$ for a 10 mg/ml solution)
Buffer	10 mM sodium phosphate, 150 mM NaCl, pH 7.2, 0.1 % $\text{NaN}_3$
Storage	2-8°C
Source	Rabbits immunized with purified rat alpha-2-macroglobulin

P.O. Box 5205 • West Chester, PA 19380 • USA