

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Rabbit Anti Mouse Ceruloplasmin (Affinity Purified)

Catalog No.	18150
Purification	Affinity purified on mouse ceruloplasmin agarose
Applications	Western blot
Protein	0.25 – 2.5 mg/ml (Assuming $A_{280} = 14$ for a 10 mg/ml solution)
Buffer	10 mM sodium phosphate, 150 mM NaCl, pH 7.2, 0.1 % NaN_3
Storage	2-8°C
Source	Serum from rabbits immunized with purified mouse ceruloplasmin

P.O. Box 5205 • West Chester, PA 19380 • USA