

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Rat Vitamin-D Binding Protein (Gc-globulin)

Catalog No.	20220
Analysis	SDS PAGE: Single ~54 kDa band
Purity	≥95%
Concentration	0.25 – 2.5 mg/ml
Buffer	150 mM NaCl, 10 mM sodium phosphate, 0.1% NaN ₃ , pH 7.2
Storage	-20 to -70°C. Avoid repeated freeze-thaws
Source	Fresh frozen rat plasma

Rat liver Vitamin-D binding protein (Gc-globulin) has recently been identified as a useful biomarker of liver injury. It is purified using a combination of affinity, gel filtration and ion-exchange chromatography.

P.O. Box 5205 • West Chester PA 19380 • USA