

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Rabbit Skeletal Muscle Troponin-C

Catalog No.	3010
Analysis	SDS PAGE, Single ~18 kDa band
Purity	≥95% by SDS PAGE
Total Protein	5 - 20 mg/ml
Buffer	150 mM sodium chloride, 10 mM sodium phosphate, 0.05% sodium azide, (pH 7.2)
Storage	-20°C to -70°C (avoid repeated freeze thawing)
Source	Rabbit skeletal muscle

Troponin-C is one of three subunits that form the troponin complex of the thin filaments of striated muscle. Skeletal muscle troponin-C is characterized by a molecular weight of ~18 kDa. It is purified using a combination of ion-exchange and affinity chromatography steps.

P.O. Box 5205 • West Chester PA 19380 • USA