

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Rabbit Cardiac Fatty Acid Binding Protein

Catalog No.	7933
Analysis	SDS PAGE: Single ~15 kDa band
Purity	≥95%
Concentration	0.2 - 2 mg/ml
Buffer	150 mM NaCl, 10 mM sodium phosphate, 0.05% NaN ₃ (pH 7.2) or 75 mM NaCl, 5 mM sodium phosphate (pH 7.2), 0.025% NaN ₃ containing 50% glycerol
Storage	-20°C (do not store at lower temperatures)
Source	Rabbit heart

Cardiac fatty acid-binding protein (FABP) is a small cytosolic protein that is abundant in the heart and skeletal muscle. Following cardiac muscle injury it rapidly appears in the blood making it a useful biomarker for detection of cardiac damage. FABP is purified using a combination of gel filtration and ion-exchange chromatography.

P.O. Box 5205 • West Chester PA 19380 • USA