

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Monkey C-Reactive Protein

Catalog No.	8103
Purity	~90% by SDS PAGE: Major ~25 kDa band under reducing conditions
Concentration	0.2 - 1 mg/ml (ELISA)
Buffer	100 mM TrisCl (pH 7.5), 200 mM NaCl, 2 mM CaCl ₂ , 0.10 % NaN ₃
Storage	-20 to -70°C. Avoid repeated freeze-thaws.
Source	Normal Rhesus monkey serum
Use Statement	For research purposes only

C-Reactive protein (CRP) is a pentameric serum protein belonging to the pentraxin family that is comprised of five identical ~25 kDa subunits. It is a positive acute phase reactant, the levels of which can increase in excess of 100-fold in response to inflammation, tissue injury or disease

P.O. Box 5205 • West Chester PA 19380 • USA