

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Rat Hemopexin***

Catalog No.	8720
Analysis	SDS PAGE: Single ~50 kDa band
Purity	≥95%
Concentration	0.5 – 2.5 mg/ml
Buffer	10 mM sodium phosphate (pH 7.2), 150 mM NaCl, 0.05% sodium azide
Storage	-20 to -70°C. Avoid repeated freeze-thaws
Source	Fresh frozen rat serum

Rat hemopexin is a serum protein of ~51 kDa. It functions as a heme scavenger in blood. It is also a positive acute phase reactant, the levels of which increase significantly as a result of inflammation, tissue injury and disease. It is purified using a combination of affinity, gel filtration and ion-exchange chromatography steps.

P.O. Box 5205 • West Chester PA 19380 • USA