

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Pig IgG

Catalog No.	30001-5
Analysis	SDS PAGE Major bands of ~50 and ~25 kDa indicative of the heavy and light chains of IgG on reducing SDS PAGE gels.
Purity	≥98%
Concentration	2 - 10 mg/ml (Assuming an extinction coefficient of 1.4 at 280 nm for a 1 mg/ml solution, path length = 1 cm)
Buffer	150 mM NaCl, 10 mM sodium phosphate, 0.1% NaN ₃ , (pH 7.2)
Storage	2 to 8°C
Source	Normal pig serum

Comments: Purified at Life Diagnostics using a combination of affinity and size exclusion chromatography

PO Box 5205 • West Chester PA 19380 • USA