

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Goat Anti-Dog C-Reactive Protein (CRP)

Catalog No.	18153
Description	Goat anti-dog CRP, affinity purified on dog-CRP agarose.
Purity	>95% by SDS PAGE
Protein	2.5-10 mg/ml (E1% at 280 nm = 14.0)
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.09% NaN ₃
Storage	2-8°C short term, ≤-20°C long term
Immunogen	Highly purified dog CRP prepared at Life Diagnostics, Inc.
Comments	Reactivity confirmed in direct ELISA using dog CRP as capture reagent and rabbit anti-goat IgG-AP for detection.

P.O. Box 5205 • West Chester, PA 19380 • USA