

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Cat C-Reactive Protein***

Catalog No.	8104
Analysis	SDS PAGE: ~20 kDa doublet on Laemmli reducing gels.
Purity	≥95%
Concentration	1 to 6 mg/ml (E <sup>1%</sup> @ 280nm = 16.53)
Buffer	50 mM Tris.Cl (pH 7.5), 200 mM NaCl, 2 mM calcium chloride, 0.05% sodium azide
Storage	-20 to -70°C. Avoid repeated freeze-thaws
Source	Pooled cat serum

Cat C-Reactive protein (CRP) is a pentameric serum protein belonging to the pentraxin family that is comprised of five ~20 kDa subunits. CRP is purified at Life Diagnostics by multiple affinity chromatography steps.

P.O. Box 5205 • West Chester PA 19380 • USA