

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Cow Serum Amyloid A***

*(serum isoform)*

Catalog No.	9007
Analysis	SDS PAGE: ~12 kDa band
Purity	90%
Concentration	2.0 mg/ml
Buffer	6M Urea, 100 mM TrisCl, 1 mM EDTA (pH 7.5)
Storage	-20 to -80°C

Serum amyloid A (SAA) is a positive acute phase reactant that is synthesized in the liver. In many species, including cows, serum levels increase significantly (>100-fold) in response to inflammation, tissue injury or disease. SAA can be diluted from the urea solution into non-urea buffers for use as assay controls and calibrators.

Rev. 091418

P.O. Box 5205 • West Chester PA 19380 • USA