

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

*Info@lifediagnostics.com • www.lifediagnostics.com*

## ***Anti-Cat Alpha-1-Acid Glycoprotein***

(Mouse Monoclonal)

Catalog No.	AGP-08-11D2
Isotype	IgG1 kappa
Analysis	SDS PAGE: Major ~60 & 30 kDa bands on reducing gels
Purity	90%
Concentration	0.25 - 2.5 mg/ml (Abs 0.1% at 280 nm = 1.40)
Buffer	10 mM sodium phosphate (pH 7.2), 150 mM NaCl, 0.09% sodium azide
Storage	2-8°C short term, -80°C long term.
Immunogen	Purified cat alpha-1-acid glycoprotein.

Comments: Developed and manufactured at Life Diagnostics, Inc. Purified from serum free media using hydroxyapatite chromatography at neutral pH. AGP-08-11D2 can be used as detection antibody with monoclonals AGP-08-4F57 or AGP-08-11B7 as capture antibodies.

**P.O. Box 5205 • West Chester PA 19380 • USA**