

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Monoclonal Anti-Dog Haptoglobin

Catalog No.	HAPT-04-9G6
Clone	9G6-1-2
Description	Mouse monoclonal IgG1k
Purity	≥90% by SDS PAGE
Protein	1 - 5 mg/ml (E1% at 280 nm = 14.0)
Applications	ELISA, Western Blot.
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.09% NaN ₃
Storage	2-8°C short term, -80°C long term
Immunogen	Dog Haptoglobin

Comments: Western blots indicate the antibody recognizes the haptoglobin heavy chain. It can be used in serial ELISA's as both capture and detection antibody. In ELISA no reactivity was observed with haptoglobin from mouse, rat, pig, cow, horse and monkey. It was developed and manufactured at Life Diagnostics, Inc. Purification was from serum free media using affinity chromatography at neutral pH.

Rev 030818

P.O. Box 5205 • West Chester, PA 19380 • USA