



Life Diagnostics, Inc.

Phone: 610-431-7707
Fax: 610-431-7818
Email: info@lifediagnostics.com
www.lifediagnostics.com

Life Diagnostics, Inc.
P.O. Box 5205
West Chester PA, 19380
USA

Monoclonal Anti-Bovine Haptoglobin

Catalog No.	HAPT-11-6A8
Description	Mouse monoclonal IgG1k
Purity	≥90% by SDS PAGE
Clone	6A8-10-26
Form	Liquid
Concentration	1 - 5 mg/ml (E ¹ % at 280 nm = 14.0)
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, 0.09% NaN ₃ , pH 7.2
Storage	2°C - 8°C short term, -80°C long term
Applications	SPARCL™, pairs with HAPT-11-3F7. Not yet tested in ELISA.
Immunogen	Purified bovine haptoglobin
Manufactured	In the United States of America

Comments: This antibody was developed using bovine haptoglobin as immunogen. It was developed and manufactured at Life Diagnostics, Inc. Purification was performed from serum-free hybridoma supernatant using molecular sieve chromatography.