

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Anti-Human LRG1 Monoclonal Clone 13H5-9***

Catalog No.	LRG1-20-13H5
Description	Anti-Human leucine-rich alpha-2-glycoprotein-1 monoclonal antibody. Purified by hydroxyapatite chromatography from serum-free culture medium. Mouse IgG1
Purity	>95% by SDS PAGE
Protein	Varies with lot.
Applications	ELISA, SPARCL. Other uses not yet investigated.
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.09% NaN <sub>3</sub>
Storage	2-8°C (up to 1 month) ≤ 70°C (longer than 1 month)
Immunogen	Human LRG1
Pairing	This antibody pairs with rabbit polyclonal anti-LRG1-20-RB

P.O. Box 5205 • West Chester, PA 19380 • USA