

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Multi-PEGylated BSA *(PEGylated with 550 Da mPEG)*

Catalog No.	PBSA-00-550
Concentration	1-5 mg/ml (Assuming an absorbance of 0.632 at 280 nm for a 1 mg/ml solution). Please note: The concentration above refers only to the BSA polypeptide concentration. It does not include the mass contributed by the 550 Da mPEG chains
Buffer	150 mM NaCl, 10 sodium phosphate, 0.09% NaN ₃ , (pH 7.2)
Storage	2-8°C short term, -20 to -80°C long term
Comments	Prepared by reaction of 550 Da mPEG succinimidyl succinate (mPEG-SS) with BSA. A 90-fold molar excess of mPEG-SS was used. Free-PEG was removed by gel-filtration.

PO Box 5205 • West Chester PA 19380 • USA