



Life Diagnostics, Inc.

Phone: 610-431-7707

Fax: 610-431-7818

Email: info@lifediagnostics.com

www.lifediagnostics.com

Life Diagnostics, Inc.

P.O. Box 5205

West Chester PA, 19380

USA

Horse Serum Amyloid A (SAA) Monoclonal

Catalog/Clone	SAA-14-10D5
Description	Mouse monoclonal anti-horse SAA Purified from serum-free culture medium Isotype: Mouse IgG1-kappa
Purity	90% by SDS PAGE
Form	Liquid
Concentration	0.5 - 2 mg/ml (Abs1% at 280 nm = 14.0)
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.1% NaN ₃
Storage	2-8°C (short term), -20°C or lower (long term)
Application	ELISA. Use as conjugate antibody paired with SAA-14-9B1 or SAA-14-16C10 as capture antibodies.
Immunogen	Horse serum amyloid A
Manufactured	USA

For research use only. This product is not designed for diagnostic procedures. This product is not for resale without authorized consent from Life Diagnostics Inc. Manufactured in the United States of America.

Revision 02/03/21