

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Cat Serum Amyloid A Monoclonal***

Catalog No.	SAA-8-13E128
Clone	13E12-8
Description	Anti-Cat SAA Purified from serum-free culture medium
Isotype	Mouse IgG2a-kappa
Purity	90% by SDS PAGE
Protein	0.5-5 mg/ml (E1% at 280 nm = 14.0)
Applications	ELISA. Use as detection antibody paired with SAA-8-15E9 as capture antibody.
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.10% NaN <sub>3</sub>
Storage	2-8°C short term, -20°C long term.
Immunogen	Cat serum amyloid A

Comments: In sandwich ELISA this works well as HRP conjugated detection antibody (~1 µg/ml) paired with SAA-8-15E9 (2.5 - 10 µg/ml).

P.O. Box 5205 • West Chester, PA 19380 • USA