

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Cat Serum Amyloid A Monoclonal

Catalog No.	SAA-8-15H9
Clone	15H9-15
Description	Anti-Cat SAA Purified from serum-free culture medium
Isotype	Mouse IgG2b-kappa
Purity	90% by SDS PAGE
Protein	0.5 - 5 mg/ml (E1% at 280 nm = 14.0)
Applications	ELISA. Use as capture antibody paired with SAA-8-13E128 as detection antibody.
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.10% NaN ₃
Storage	2-8°C short term, -20°C long term.
Immunogen	Cat serum amyloid A

Comments: In sandwich ELISA this works well as coating antibody (2.5 – 10 µg/ml) paired with HRP-conjugated SAA-8-13E128 (~ 1 µg/ml).

P.O. Box 5205 • West Chester, PA 19380 • USA