



**Veterinary
Biomarkers, Inc.**

Phone: 484-200-9651
Email: info@vetbiomarkers.com
www.vetbiomarkers.com

Veterinary Biomarkers, Inc.
P.O. Box 2532
West Chester PA, 19380
USA

Cow SAA3 Monoclonal Milk SAA Isoform Clone 5H6-44-12

Catalog No.	SAA3-11-5H6
Lot	D-E2219A
Description	Mouse anti-Cow SAA3 Purified from serum-free culture medium by molecular sieve and hydroxyapatite chromatography
Purity	≥95% by SDS PAGE
Form	Liquid
Concentration	0.5 mg/ml (E ^{1%} at 280 nm = 14.0)
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, 0.09% NaN ₃ , pH 7.2
Storage	2°C - 8°C
Recertification	May 2029
Immunogen	SAA3 Peptide
Applications	ELISA and SPARCL™
Comments	Pairs with SAA3-11-11H3

For research use only. This product is not designed for diagnostic procedures. This product is not for resale without authorized consent from Life Diagnostics Inc. Manufactured in the United States of America.

Revision 09172019