



Life Diagnostics, Inc.

Phone: 610-431-7707

Fax: 610-431-7818

Email: [info@lifediagnostics.com](mailto:info@lifediagnostics.com)

[www.lifediagnostics.com](http://www.lifediagnostics.com)

Life Diagnostics, Inc.

P.O. Box 5205

West Chester PA, 19380

USA

## Stabilized Human Cardiac Troponin-I

Catalog No.	1230-R
Form	Liquid
Concentration	40 µg/ml Troponin-I
Buffer	500 mM NaCl, 5 mM sodium phosphate, 2 mM CaCl <sub>2</sub> , 0.5 mM EDTA 0.05% NaN <sub>3</sub> , pH 7.0 containing recombinant rabbit skeletal muscle troponin-C as a stabilizer
Storage	-70°C Avoid repeated freeze-thaws
Source	Recombinant Human Cardiac Troponin-I Recombinant Rabbit SkM Troponin-C
Comments	Prepare stock dilutions in delipidated human serum containing 1 mM CaCl <sub>2</sub> . Freeze as aliquots and store at -70°C. Stability in serum varies with batch.
Manufactured	USA