

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Pig Cardiac Fatty Acid Binding Protein***

Catalog No.	7935
Analysis	SDS PAGE: Single ~15 kDa band
Purity	≥95%
Concentration	0.2 - 2 mg/ml
Buffer	150 mM NaCl, 10 mM sodium phosphate, 0.05% NaN <sub>3</sub> (pH 7.2)
Storage	-20°C
Source	Pig heart

Cardiac fatty acid-binding protein (FABP) is a small cytosolic protein that is abundant in the heart and skeletal muscle. Following cardiac muscle injury it rapidly appears in the blood making it a useful biomarker for detection of cardiac damage. FABP is purified using a combination of gel filtration and ion-exchange chromatography.

P.O. Box 5205 • West Chester PA 19380 • USA