

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Monkey IgM

Catalog No.	30000-4
Analysis	Major bands of ~70 and ~25 kDa indicative of the heavy and light chains of IgM on reducing SDS PAGE gels. IgG and IgA levels \leq 2%
Purity	\geq 90%
Concentration	0.5 - 2 mg/ml (Assuming an extinction coefficient of 1.4 at 280 nm for a 1 mg/ml solution, path length = 1 cm)
Buffer	150 mM NaCl, 10 mM sodium phosphate, 0.1% NaN ₃ , (pH 7.2)
Storage	2 to 8°C
Source	Normal monkey serum (Rhesus)

Comments: Purified at Life Diagnostics using size exclusion chromatography

PO Box 5205 • West Chester PA 19380 • USA