

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Rat IgA***

Catalog No.	30002-2
Analysis	Major bands of ~60 and ~25 kDa indicative of the heavy and light chains of IgA on reducing SDS PAGE gels.
Purity	90%
Concentration	0.5 - 2 mg/ml (Assuming an extinction coefficient of 1.4 at 280 nm for a 1 mg/ml solution, path length = 1 cm)
Buffer	150 mM NaCl, 10 mM sodium phosphate, 0.1% NaN <sub>3</sub> , (pH 7.2)
Storage	2 to 8°C
Source	Normal rat serum

Comments: Purified at Life Diagnostics using a combination of affinity and size exclusion chromatography

PO Box 5205 • West Chester PA 19380 • USA