

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Monkey IgA

| | |
|---------------|---|
| Catalog No. | 30002-4 |
| Analysis | Major bands of ~60 and ~25 kDa indicative of the heavy and light chains of IgA on reducing SDS PAGE gels. IgG levels \leq 2.5% by ELISA |
| Purity | \geq 90% |
| Concentration | 0.5 - 2 mg/ml (Assuming an extinction coefficient of 1.4 at 280 nm for a 1 mg/ml solution, path length = 1 cm) |
| Buffer | 150 mM NaCl, 10 mM sodium phosphate, 0.1% NaN ₃ , (pH 7.2) |
| Storage | 2 to 8°C |
| Source | Normal monkey serum (Rhesus) |

Comments: Purified at Life Diagnostics using a combination of affinity and size exclusion chromatography

PO Box 5205 • West Chester PA 19380 • USA