

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Mono PEGylated BSA***

***(PEGylated with 20kDa pentanoic acid-PEG maleimide)***

Catalog No.	PBSA-01P
Analysis	SDS PAGE & HPLC
Purity	≥90%
Concentration	1 - 10 mg/ml
Buffer	150 mM NaCl, 10 sodium phosphate, 0.1% NaN <sub>3</sub> , (pH 7.2)
Storage	-80°C
Comments	Prepared by conjugation of pentanoic acid-PEG-maleimide (20 kDa) to native BSA followed by gel-filtration chromatography to remove free PEG and unconjugated BSA.

PO Box 5205 • West Chester PA 19380 • USA