

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

PEG-BSA Coated Plate

(12 x 8-well strips)

Catalog No.	PBSA20PL
Analysis	ELISA. CV's \leq 10%
Coating	Purified mono mPEGylated BSA (20 kDa mPEG)
Blocking	Blocked and stabilized using a proprietary PEG-free buffer (PEGB125-1).
Storage	-20°C (maintain desiccated)

Comments: For use in the detection and quantitation of anti-PEG antibodies. All buffers should be free of extraneous PEGs and PEGylated compounds including detergents such as Tween-20. We recommend use of our PEG-free dilution (PEGD50-1) and wash (PEGW50-20) buffers.

Rev 062013

P.O. Box 5205 • West Chester PA 19380 • USA