

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

Sheep C-Reactive Protein

Catalog No.	8107
Analysis	SDS PAGE: ~25 kDa band on Laemmli reducing gels.
Purity	≥90%
Concentration	1-10 mg/ml ($E^{1\%}$ @ 280nm = 19.57)
Buffer	10 mM Tris.Cl (pH 8.0), 150 mM NaCl, 1 mM EDTA, 0.09% sodium azide
Storage	-20 to -70°C. Avoid repeated freeze-thaws
Source	Fresh frozen sheep serum

Sheep C-Reactive protein (CRP) is a pentameric serum protein belonging to the pentraxin family that is comprised of five ~22 kDa subunits. It is purified at Life Diagnostics by affinity chromatography.

P.O. Box 5205 • West Chester PA 19380 • USA