

Life Diagnostics, Inc.

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

PEG Dilution Buffer

Catalog No.	PEGD50-1
Analysis	Not applicable
Composition	Proprietary PEG-free buffer containing 0.05% ProClin 300 as preservative.
Volume	50 ml
Storage	2-8°C

Comments: Dilution buffer for use in PEG ELISA's and western blots, This product can be used together with our PEG-BSA coated plates (PBSA20PL), PEG Blocking buffer (PEGB125-1) and PEG wash buffer (PEGW50-20).

Please Note: This is not intended for use with the MPEG ELISA kit

Rev 082015

P.O. Box 5205 • West Chester PA 19380 • USA