

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## **PEG Wash Buffer, 20x**

(50 ml)

Catalog No.	PEGW50-20
Analysis	Not applicable
Composition	Proprietary PEG-free buffer containing 0.05% ProClin 300 as preservative.
Instructions	Dilute with 19-volumes of distilled or deionized water prior to use.
Storage	2-8°C

Comments: Wash buffer concentrate (20 x) for use in ELISA or western-blot detection of PEG, PEGylated proteins and anti-PEG antibodies. This product can be used together with our PEG-BSA coated plates (PBSA20PL) and PEG-free dilution buffer (PEGD50-1).

Please Note: This is not intended for use with the MPEG ELISA kit

Rev 082015

P.O. Box 5205 • West Chester PA 19380 • USA