



Phone: 610-431-7707  
Fax: 610-431-7818  
Email: [info@lifediagnostics.com](mailto:info@lifediagnostics.com)  
[www.lifediagnostics.com](http://www.lifediagnostics.com)

Life Diagnostics, Inc.  
P.O. Box 5205  
West Chester PA, 19380  
USA

## Horse Serum Amyloid A (SAA) Monoclonal

Catalog/Clone	SAA-14-9B1
Description	Mouse monoclonal anti-horse SAA Purified from serum-free culture medium Isotype: Mouse IgG2b-kappa
Purity	90% by SDS PAGE
Form	Liquid
Concentration	0.5 - 2 mg/ml (Abs1% at 280 nm = 14.0)
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.2, 0.1% NaN <sub>3</sub>
Storage	2-8°C (short term), -20°C or lower (long term)
Application	ELISA. Use as capture antibody paired with SAA-14-10D5 conjugate antibody.
Immunogen	Horse serum amyloid A
Manufactured	USA

For research use only. This product is not designed for diagnostic procedures. This product is not for resale without authorized consent from Life Diagnostics Inc. Manufactured in the United States of America.

Revision 02/03/21